

ARCONIC - Aluminium base 6063, Class: RM, Alloy: 6063, cast, typical analysis, ppm values not shown, Disc 64 mm x 25 mm

Art. ID ALCOA-WA-6063
Unit disc
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Copper (Cu)	[7440-50-8]	0,0105	%		
	Iron (Fe)	[7439-89-6]	0,164	%		
	Gallium (Ga)	[7440-55-3]	0,014	%		
	Magnesium (Mg)	[7439-95-4]	0,482	%		
	Manganese (Mn)	[7439-96-5]	0,0315	%		
	Silicon (Si)	[7440-21-3]	0,45	%		
	Titanium (Ti)	[7440-32-6]	0,0133	%		
	Vanadium (V)	[7440-62-2]	0,01	%		