

### **Custom Solution (3 components) in Acetone**

Art. ID LM24-CUS-70093  
Unit 4,5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1090 II E2

#### Description

Stored under normal laboratory conditions. Conservation 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Phthalic acid, bis-phenyl ester (Bis-phenyl phthalate)	[84-62-8]	100	µg/mL		
	Di-n-octyl phthalate	[117-84-0]	100	µg/mL		
	Phthalic acid, bis-decyl ester (Bis-decyl phthalate)	[84-77-5]	100	µg/mL		