

## **High-purity zinc, Block**

Art. ID	NIST-682
Unit	disc
Deliverydetails	No Dangerous Good /not restricted

### Description

This Standard Reference Material (SRM®) is intended primarily for the calibration of instruments and the evaluation of chemical methods used in zinc analysis. NIST-682 is in the form of a semicircular bar segment, 57 mm diameter (2 1/4 inch), 25,4 mm (1 inch) deep at mid-diameter and 19 mm long (3/4 inch). /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Copper (Cu)	[7440-50-8]	0,0000042	%		