

3-Amino-2,7-naphthalenedisulfonic Acid Sodium Salt

Art. ID TRC-A610265-25G
Unit 25 g
Deliverydetails to be checked

Description

Application Note: 3-Amino-2,7-naphthalenedisulfonic Acid Sodium Salt, can be used in various reactions including in the synthesis of 2,3,6-Trichloronaphthalene (T774070), that binds to the aryl hydrocarbon receptor (AhR) and induce a wide range of pleiotrophic effects.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	3-Amino-2,7-naphthalene disulfonic Acid Sodium Salt	[5332-41-2]				