

(4-Acetyl-piperazin-1-yl)acetic Acid

Art. ID TRC-A187238-250MG
Unit 250 mg

Description

Application Note: (4-Acetyl-Piperazin-1-Yl)-Acetic Acid is a useful research chemical. It is used in the preparation of phenylpiperidines as neurokinin NK1 and NK3 antagonists for treatment of neurokinin-mediated diseases

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	(4-Acetyl-piperazin-1-y l)acetic Acid	[705941-45-3]				