

L-Arginine-L-phenylalanine-L-arginine-L-tryptophan-L-proline-L-tyrosine-L-arginine-L-isoleucine-L-arginine-L-glutamic acid-L-phenylalanine TFA Salt (95%)

Art. ID TRC-A169125-25MG
Unit 25 mg
Deliverydetails to be checked

Description

Application Note:

L-Arginine-L-phenylalanine-L-arginine-L-tryptophan-L-proline-L-tyrosine-L-arginine-L-isoleucine-L-arginine-L-glutamic acid-L-phenylalanine is a useful peptide compound.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	L-Arginine-L-phenylalanine-L-arginine-L-tryptophan-L-proline-L-tyrosine-L-arginine-L-isoleucine-L-arginine-L-glutamic acid-L-phenylalanine TFA Salt					