

2-Isobutyl-3-methoxypyrazine 100 µg/mL in Methanol, analytical standard (Supelco)

Art. ID SAF-47528-U
Unit 1 x 1 mL
Deliverydetails De Minimis UN1230 (6.1) II E2

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|----------------------------------|--------------|---------------------|-------|--------|--------|
| | 2-Isobutyl-3-methoxypyr azine | [24683-00-9] | 100 | µg/mL | | |