

### **Custom Solution - Aluminium Al - 10 g/L in diluted HNO3 for ICP**

Art. ID LM24-CUS-70600  
Unit 125 mL bottle  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

#### Description

Conservation: 24 months (originally closed bottle)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	10	g/L		