

Cast Iron, Ni-Hard, Type IV

Art. ID	NIST-892
Unit	150 g
Deliverydetails	No Dangerous Good /not restricted

Description

This Standard Reference Material (SRM®) is in the form of small ganules. Certified values /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Carbon (C)	[7440-44-0]	3,33	%		
	Manganese (Mn)	[7439-96-5]	0,76	%		
	Phosphorus (P)	[7723-14-0]	0,054	%		
	Sulfur (S)	[7704-34-9]	0,015	%		
	Silicon (Si)	[7440-21-3]	1,83	%		
	Copper (Cu)	[7440-50-8]	0,207	%		
	Nickel (Ni)	[7440-02-0]	5,53	%		
	Chromium (Cr)	[7440-47-3]	10,18	%		
	Vanadium (V)	[7440-62-2]	0,041	%		
	Molybdenum (Mo)	[7439-98-7]	0,2	%		
	Cobalt (Co)	[7440-48-4]	0,31	%		