

**Aluminium base, Cu/Si/Zn, RM, cast, Mushroom 58 mm x 7 mm**

Art. ID IMI-24WM4  
Unit each  
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Chromium (Cr)     | [7440-47-3] | 0,045               | %    |        |        |
|                  | Copper (Cu)       | [7440-50-8] | 2,93                | %    |        |        |
|                  | Iron (Fe)         | [7439-89-6] | 0,69                | %    |        |        |
|                  | Magnesium (Mg)    | [7439-95-4] | 0,26                | %    |        |        |
|                  | Manganese (Mn)    | [7439-96-5] | 0,43                | %    |        |        |
|                  | Nickel (Ni)       | [7440-02-0] | 0,45                | %    |        |        |
|                  | Lead (Pb)         | [7439-92-1] | 0,24                | %    |        |        |
|                  | Silicon (Si)      | [7440-21-3] | 8,85                | %    |        |        |
|                  | Tin (Sn)          | [7440-31-5] | 0,16                | %    |        |        |
|                  | Titanium (Ti)     | [7440-32-6] | 0,16                | %    |        |        |
|                  | Zinc (Zn)         | [7440-66-6] | 1,5                 | %    |        |        |