

Aluminium base, Cu/Fe/Si, RM, cast, Mushroom 58 mm x 7 mm

Art. ID IMI-24WM6
Unit each
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Chromium (Cr)	[7440-47-3]	0,19	%		
	Copper (Cu)	[7440-50-8]	1,85	%		
	Iron (Fe)	[7439-89-6]	1,02	%		
	Magnesium (Mg)	[7439-95-4]	0,2	%		
	Manganese (Mn)	[7439-96-5]	0,63	%		
	Nickel (Ni)	[7440-02-0]	0,45	%		
	Lead (Pb)	[7439-92-1]	0,26	%		
	Silicon (Si)	[7440-21-3]	6,55	%		
	Tin (Sn)	[7440-31-5]	0,1	%		
	Titanium (Ti)	[7440-32-6]	0,23	%		
	Zinc (Zn)	[7440-66-6]	0,56	%		