

Aluminium base, Cu/Fe/Si, RM, cast, Mushroom 58 mm x 7 mm

Art. ID IMI-24WM6
Unit each
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Chromium (Cr) | [7440-47-3] | 0,19 | % | | |
| | Copper (Cu) | [7440-50-8] | 1,85 | % | | |
| | Iron (Fe) | [7439-89-6] | 1,02 | % | | |
| | Magnesium (Mg) | [7439-95-4] | 0,2 | % | | |
| | Manganese (Mn) | [7439-96-5] | 0,63 | % | | |
| | Nickel (Ni) | [7440-02-0] | 0,45 | % | | |
| | Lead (Pb) | [7439-92-1] | 0,26 | % | | |
| | Silicon (Si) | [7440-21-3] | 6,55 | % | | |
| | Tin (Sn) | [7440-31-5] | 0,1 | % | | |
| | Titanium (Ti) | [7440-32-6] | 0,23 | % | | |
| | Zinc (Zn) | [7440-66-6] | 0,56 | % | | |