

**Aluminium base, Cu/Si/Zn, RM, cast, Mushroom 58 mm x 7 mm**

Art. ID IMI-2WM2  
Unit each  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Chromium (Cr)	[7440-47-3]	0,05	%		
	Copper (Cu)	[7440-50-8]	1,1	%		
	Iron (Fe)	[7439-89-6]	0,85	%		
	Magnesium (Mg)	[7439-95-4]	0,095	%		
	Manganese (Mn)	[7439-96-5]	0,48	%		
	Nickel (Ni)	[7440-02-0]	0,3	%		
	Lead (Pb)	[7439-92-1]	0,15	%		
	Silicon (Si)	[7440-21-3]	9,76	%		
	Tin (Sn)	[7440-31-5]	0,14	%		
	Titanium (Ti)	[7440-32-6]	0,16	%		
	Zinc (Zn)	[7440-66-6]	1,2	%		