

**Aluminium base, Cu/Fe/Si, RM, cast, Mushroom 58 mm x 7 mm**

Art. ID IMI-4P5  
Unit each  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Chromium (Cr)	[7440-47-3]	0,19	%		
	Copper (Cu)	[7440-50-8]	4,25	%		
	Iron (Fe)	[7439-89-6]	1,03	%		
	Magnesium (Mg)	[7439-95-4]	<0,005	%		
	Manganese (Mn)	[7439-96-5]	0,14	%		
	Nickel (Ni)	[7440-02-0]	<0,005	%		
	Lead (Pb)	[7439-92-1]	<0,005	%		
	Silicon (Si)	[7440-21-3]	3,94	%		
	Tin (Sn)	[7440-31-5]	<0,005	%		
	Titanium (Ti)	[7440-32-6]	<0,005	%		
	Zinc (Zn)	[7440-66-6]	0,61	%		