

### **4-Amino-2-methyl-5-pyrimidinethanesulfonic acid**

Art. ID TRC-A578025-250MG

Unit 250 mg

Deliverydetails to be checked

#### Description

Application Note: 4-Amino-2-methyl-5-pyrimidinethanesulfonic acid is used in the biosynthesis of Thiamine (T344185) in yeast.

| Text/Information | Analyte/Parameter                                       | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|---------------------------------------------------------|-------------|---------------------|------|--------|--------|
|                  | 4-Amino-2-methyl-5-pyri<br>midinethanesulfonic a<br>cid | [2908-73-8] |                     |      |        |        |