

N-[[[(1R)-2-(6-Amino-9H-purin-9-yl)-1-methylethoxy]methyl]hydroxyphosphinyl]-D-alanine

Art. ID TRC-A612560-100MG

Unit 100 mg

Deliverydetails to be checked

Description

Application Note: N-[[[(1R)-2-(6-Amino-9H-purin-9-yl)-1-methylethoxy]methyl]hydroxyphosphinyl]-D-alanine is a metabolite of Tenofovir (T018500), an acyclic phosphonate nucleotide analogue and reverse transcriptase inhibitor. It is used as an anti-HIV agent. Antiviral.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	N-[[[(1R)-2-(6-Amino-9H -purin-9-yl)-1-methylet hoxy]methyl]hydroxyphos phinyl]-D-alanine	[376633-24-8]				