

(3S)-3-Amino-5-hexynoic acid

Art. ID TRC-A612570-2.5G
Unit 2,5 g
Deliverydetails to be checked

Description

Application Note: (3S)-3-Amino-5-hexynoic acid is an unusual beta-3-homo-amino acid and can be analyzed using direct and indirect reversed-phase (RP) high-performance liquid chromatography (HPLC).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	(3S)-3-Amino-5-hexynoic acid	[270596-46-8]				