

Aluminium base, RM, 5001, wrought, Disc 57 mm x 25 mm

Art. ID ALCAN-5001-BI
Unit disc
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Beryllium (Be) | [7440-41-7] | 0,0012 | % | | |
| | Bismuth (Bi) | [7440-69-9] | 0,052 | % | | |
| | Calcium (Ca) | [7440-70-2] | 0,0018 | % | | |
| | Cadmium (Cd) | [7440-43-9] | 0,0013 | % | | |
| | Chromium (Cr) | [7440-47-3] | 0,049 | % | | |
| | Copper (Cu) | [7440-50-8] | 19,5 | % | | |
| | Iron (Fe) | [7439-89-6] | 0,3 | % | | |
| | Gallium (Ga) | [7440-55-3] | 0,02 | % | | |
| | Magnesium (Mg) | [7439-95-4] | 0,053 | % | | |
| | Manganese (Mn) | [7439-96-5] | 0,05 | % | | |
| | Nickel (Ni) | [7440-02-0] | 0,052 | % | | |
| | Lead (Pb) | [7439-92-1] | 0,048 | % | | |
| | Antimony (Sb) | [7440-36-0] | 0,004 | % | | |
| | Silicon (Si) | [7440-21-3] | 0,36 | % | | |
| | Tin (Sn) | [7440-31-5] | 0,047 | % | | |
| | Strontium (Sr) | [7440-24-6] | 0,0027i | % | | |
| | Titanium (Ti) | [7440-32-6] | 0,049 | % | | |
| | Vanadium (V) | [7440-62-2] | 0,025 | % | | |
| | Zinc (Zn) | [7440-66-6] | 0,045 | % | | |
| | Zirconium (Zr) | [7440-67-7] | 0,003 | % | | |