

**Phenol D6 2000 µg/mL in Methylene chloride (Dichloromethane), analytical standard (Supelco)**

Art. ID                    SAF-48716-U  
Unit                        1 x 1 mL  
Deliverydetails        De Minimis UN1593 III E1

| Text/Information | Analyte/Parameter | CAS number   | Concentration/Value | Unit  | Method | Source |
|------------------|-------------------|--------------|---------------------|-------|--------|--------|
|                  | Phenol D6         | [13127-88-3] | 2000                | µg/mL |        |        |