

**Aluminium base, RM, A333, cast, Disc 57 mm x 25 mm**

Art. ID                   ALCAN-A333.0-AF  
Unit                     disc  
Deliverydetails        No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Beryllium (Be)	[7440-41-7]	<0,0001	%		
	Bismuth (Bi)	[7440-69-9]	0,021	%		
	Calcium (Ca)	[7440-70-2]	0,0043i	%		
	Chromium (Cr)	[7440-47-3]	0,015	%		
	Copper (Cu)	[7440-50-8]	3,54	%		
	Iron (Fe)	[7439-89-6]	0,85	%		
	Gallium (Ga)	[7440-55-3]	0,03	%		
	Lithium (Li)	[7439-93-2]	0,0015i	%		
	Magnesium (Mg)	[7439-95-4]	0,082i	%		
	Manganese (Mn)	[7439-96-5]	0,35	%		
	Sodium (Na)	[7440-23-5]	0,0033i	%		
	Nickel (Ni)	[7440-02-0]	0,13	%		
	Phosphorus (P)	[7723-14-0]	0,001	%		
	Lead (Pb)	[7439-92-1]	0,015	%		
	Antimony (Sb)	[7440-36-0]	0,0057	%		
	Silicon (Si)	[7440-21-3]	9,13	%		
	Tin (Sn)	[7440-31-5]	0,017	%		
	Strontium (Sr)	[7440-24-6]	0,0029i	%		
	Titanium (Ti)	[7440-32-6]	0,032	%		
	Vanadium (V)	[7440-62-2]	0,03	%		
	Zinc (Zn)	[7440-66-6]	2,34	%		
	Zirconium (Zr)	[7440-67-7]	0,018	%		