

Aluminium base, RM, A333+Cr, cast, Disc 57 mm x 25 mm

Art. ID ALCAN-A333.0-AI
Unit disc
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Beryllium (Be)	[7440-41-7]	0,0012	%		
	Bismuth (Bi)	[7440-69-9]	0,019	%		
	Calcium (Ca)	[7440-70-2]	0,0021	%		
	Cadmium (Cd)	[7440-43-9]	0,0012	%		
	Chromium (Cr)	[7440-47-3]	0,11	%		
	Copper (Cu)	[7440-50-8]	3,47	%		
	Iron (Fe)	[7439-89-6]	0,9	%		
	Gallium (Ga)	[7440-55-3]	0,017	%		
	Lithium (Li)	[7439-93-2]	0,0022i	%		
	Magnesium (Mg)	[7439-95-4]	0,14	%		
	Manganese (Mn)	[7439-96-5]	0,33	%		
	Sodium (Na)	[7440-23-5]	0,0028i	%		
	Nickel (Ni)	[7440-02-0]	0,15	%		
	Phosphorus (P)	[7723-14-0]	0,0032i	%		
	Lead (Pb)	[7439-92-1]	0,13	%		
	Antimony (Sb)	[7440-36-0]	0,0039	%		
	Silicon (Si)	[7440-21-3]	9,7	%		
	Tin (Sn)	[7440-31-5]	0,17	%		
	Strontium (Sr)	[7440-24-6]	0,0031i	%		
	Titanium (Ti)	[7440-32-6]	0,12	%		
	Vanadium (V)	[7440-62-2]	0,02	%		
	Zinc (Zn)	[7440-66-6]	2,56	%		
	Zirconium (Zr)	[7440-67-7]	0,0029	%		