

Aluminium base, RM, 3003 M, wrought, Disc 57 mm x 25 mm

Art. ID ALCAN-X3003-AD
Unit disc
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Beryllium (Be)	[7440-41-7]	0,0011	%		
	Bismuth (Bi)	[7440-69-9]	0,0064	%		
	Cadmium (Cd)	[7440-43-9]	0,0011	%		
	Chromium (Cr)	[7440-47-3]	0,015	%		
	Copper (Cu)	[7440-50-8]	0,16	%		
	Iron (Fe)	[7439-89-6]	0,53	%		
	Gallium (Ga)	[7440-55-3]	0,021	%		
	Magnesium (Mg)	[7439-95-4]	0,014	%		
	Manganese (Mn)	[7439-96-5]	1,17	%		
	Nickel (Ni)	[7440-02-0]	0,018	%		
	Phosphorus (P)	[7723-14-0]	0,003	%		
	Lead (Pb)	[7439-92-1]	0,0077	%		
	Antimony (Sb)	[7440-36-0]	0,0025	%		
	Silicon (Si)	[7440-21-3]	0,37	%		
	Tin (Sn)	[7440-31-5]	0,018	%		
	Titanium (Ti)	[7440-32-6]	0,021	%		
	Vanadium (V)	[7440-62-2]	0,019	%		
	Zinc (Zn)	[7440-66-6]	1,39	%		
	Zirconium (Zr)	[7440-67-7]	0,0024	%		