

## **Out of Stock - Item is not available at this time - Creatinine in Frozen Human Serum**

Art. ID                      NIST-967a  
Unit                        set(4) (2 each conc.)  
Deliverydetails           Dry ice shipment /not restricted

### Description

This Standard Reference Material (SRM®) is intended primarily for use in evaluating the accuracy of procedures for the determination of creatinine in human serum. It is also intended for use in validating working or secondary reference materials. A unit of NIST-967a consists of four stoppered vials of frozen human serum, two vials each at two different creatinine concentration levels. Each vial contains 1.0 mL of human serum. /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Level 1	Creatinine	[60-27-5]	0,0749	mmol/L		
Level 1	Creatinine	[60-27-5]	0,847	mg/dL		
Level 2	Creatinine	[60-27-5]	0,3427	mmol/L		
Level 2	Creatinine	[60-27-5]	3,877	mg/dL		