

**Out of Stock - Item is not available at this time - Boric acid - Acidimetric standard**

Art. ID                      NIST-973  
Unit                         100 g  
Deliverydetails            No Dangerous Good /not restricted

Description

This Standard Reference Material (SRM®) is certified as a chemical of known assay and is intended for use as a primary acidimetric standard. The SRM® consists of highly purified boric acid (H<sub>3</sub>BO<sub>3</sub>) in a 100 g bottle. /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	wC <sub>6</sub> H <sub>5</sub> COOH		100.009 % ± 0.010 %			