

Aluminium base, CRM, 5251, cast, Disc 55 x 30 mm

Art. ID S-PY-6036
Unit disc
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silicon (Si)	[7440-21-3]	0,102 ± 0,003	%		
	Iron (Fe)	[7439-89-6]	0,276 ± 0,009	%		
	Copper (Cu)	[7440-50-8]	0,0065 ± 0,0004	%		
	Manganese (Mn)	[7439-96-5]	0,383 ± 0,006	%		
	Magnesium (Mg)	[7439-95-4]	2,07 ± 0,04	%		
	Chromium (Cr)	[7440-47-3]	0,1248 ± 0,0031	%		
	Nickel (Ni)	[7440-02-0]	0,0487 ± 0,0011	%		
	Zinc (Zn)	[7440-66-6]	0,0936 ± 0,0013	%		
	Titanium (Ti)	[7440-32-6]	0,0527 ± 0,0019	%		
	Beryllium (Be)	[7440-41-7]	0,00096 ± 0,00007	%		
	Sodium (Na)	[7440-23-5]	0,00034 ± 0,00005	%		
	Lead (Pb)	[7439-92-1]	0,0195 ± 0,0009	%		
	Tin (Sn)	[7440-31-5]	0,0023 ± 0,0004	%		
	Zirconium (Zr)	[7440-67-7]	0,0504 ± 0,0011	%		