

3-Acetyl-4-hydroxy-7,8-dimethyl-2(1H)-quinolinone

Art. ID TRC-A178900-1G

Unit 1 g

Deliverydetails to be checked

Description

Application Note: 3-Acetyl-4-hydroxy-7,8-dimethyl-2(1H)-quinolinone is an intermediate used in the synthesis of Repirinast-d4 (R144562), which is labelled Repirinast (R144560). Antiallergic.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	3-Acetyl-4-hydroxy-7,8-dimethyl-2(1H)-quinolinone	[63768-46-7]				