

## **2-Amino-acetonitrile-2,2-d2 Hydrochloride**

Art. ID TRC-A580057-250MG  
Unit 250 mg  
Deliverydetails to be checked

### Description

Application Note: Isotope labelled Aminoacetonitrile Hydrochloride (A580055), a useful synthetic intermediate. It is a reagent used to synthesize dipeptide nitriles as reversible and potent cathepsin S inhibitors. It can also be used to prepare substituted cyclic ureas as possible HIV protease inhibitors.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2-Amino-acetonitrile-2, 2-d2 Hydrochloride	[1221897-87-5]				