

Steel Charpy V-notch test pieces (ISO 148-1)

| | |
|-----------------|-----------------------------------|
| Art. ID | ERM-FA013bk |
| Unit | set of 5 bars |
| Deliverydetails | No Dangerous Good /not restricted |

Description

A unit consists of five Charpy V-notch test pieces, which are rectangular steel bars of nominal dimensions 55 mm x 10 mm x 10 mm, with one V-notch, accurately machined to tolerances imposed in ISO 148-3

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|----------------------|------------|---------------------|------|--------|--------|
| at 20 ± 2 °C | Absorbed energy (KV) | | 27,8 | J | | |
| at 0 ± 2°C | Absorbed energy (KV) | | 26,1 | J | | |