

## **Certified Reference Material - Total Cyanide QC plus QC Standard**

Art. ID                      ERA-4093  
Unit                        each  
Deliverydetails            ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1824 III E1

### Description

One screw-cap vial yields up to 2 liters after dilution. /// The certified values vary from batch to batch /// Methods: EPA 335.2, EPA 335.3, EPA 9010B, EPA 9012A, SM 18/19thED 4500 CN-F, SM 18/19thED 4500CN-E, SM 18th ED 4500-CN D, SM 18th ED 4500-CN E, SM 18th ED 4500-CN F

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Total cyanide		1,0 -5,0	mg/L		