

Aluminium base, Si, CRM, cast, Disc 50 x 30 mm

Art. ID S-PY-9331
Unit disc
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silicon (Si)	[7440-21-3]	9,68 ± 0,17	%		
	Iron (Fe)	[7439-89-6]	0,129 ± 0,002	%		
	Copper (Cu)	[7440-50-8]	0,0021 ± 0,0001	%		
	Manganese (Mn)	[7439-96-5]	0,925 ± 0,012	%		
	Magnesium (Mg)	[7439-95-4]	0,6082 ± 0,0142	%		
	Chromium (Cr)	[7440-47-3]	0,0007 ± 0,0001	%		
	Nickel (Ni)	[7440-02-0]	0,0013 ± 0,0001	%		
	Zinc (Zn)	[7440-66-6]	0,0048 ± 0,0002	%		
	Titanium (Ti)	[7440-32-6]	0,0122 ± 0,0004	%		
	Lead (Pb)	[7439-92-1]	0,0194 ± 0,0006	%		
	Antimony (Sb)	[7440-36-0]	0,0942 ± 0,0024	%		
	Tin (Sn)	[7440-31-5]	0,0186 ± 0,0025	%		