

Aluminium base, Si, SUS, AlSi7Mg Monitor, cast, Disc 50 x 50 mm

Art. ID S-PY-9520
Unit disc
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silicon (Si)	[7440-21-3]	6,6	%		
	Iron (Fe)	[7439-89-6]	0,1	%		
	Copper (Cu)	[7440-50-8]	0,012	%		
	Manganese (Mn)	[7439-96-5]	0,005	%		
	Magnesium (Mg)	[7439-95-4]	0,34	%		
	Nickel (Ni)	[7440-02-0]	0,003	%		
	Zinc (Zn)	[7440-66-6]	0,017	%		
	Titanium (Ti)	[7440-32-6]	0,12	%		
	Calcium (Ca)	[7440-70-2]	0,0044	%		
	Sodium (Na)	[7440-23-5]	0,0005	%		
	Phosphorus (P)	[7723-14-0]	0,0009	%		
	Lead (Pb)	[7439-92-1]	<0,0001	%		
	Antimony (Sb)	[7440-36-0]	<0,008	%		
	Tin (Sn)	[7440-31-5]	0,0004	%		
	Strontium (Sr)	[7440-24-6]	0,052	%		