

**Aluminium base, Mn, SUS, 3003 Monitor, cast, Disc 50 x 50 mm**

Art. ID S-PY-9627  
Unit disc  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silicon (Si)	[7440-21-3]	0,2	%		
	Iron (Fe)	[7439-89-6]	0,57	%		
	Copper (Cu)	[7440-50-8]	0,13	%		
	Manganese (Mn)	[7439-96-5]	1,06	%		
	Magnesium (Mg)	[7439-95-4]	0,0004	%		
	Chromium (Cr)	[7440-47-3]	0,0223	%		
	Nickel (Ni)	[7440-02-0]	0,001	%		
	Zinc (Zn)	[7440-66-6]	0,057	%		
	Titanium (Ti)	[7440-32-6]	0,022	%		
	Cadmium (Cd)	[7440-43-9]	<0,0001	%		
	Mercury (Hg)	[7439-97-6]	0,0001	%		
	Lithium (Li)	[7439-93-2]	0,00002	%		
	Sodium (Na)	[7440-23-5]	<0,00002	%		
	Zirconium (Zr)	[7440-67-7]	0,01	%		