

Aluminium base, ROHS/WEEE, SUS, 1100 Monitor, cast, Disc 50 x 50 mm

Art. ID S-PY-9630
Unit disc
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silicon (Si)	[7440-21-3]	0,1	%		
	Iron (Fe)	[7439-89-6]	0,46	%		
	Copper (Cu)	[7440-50-8]	0,062	%		
	Manganese (Mn)	[7439-96-5]	0,0123	%		
	Magnesium (Mg)	[7439-95-4]	0,0006	%		
	Chromium (Cr)	[7440-47-3]	0,0216	%		
	Nickel (Ni)	[7440-02-0]	0,008	%		
	Zinc (Zn)	[7440-66-6]	0,054	%		
	Titanium (Ti)	[7440-32-6]	0,018	%		
	Cadmium (Cd)	[7440-43-9]	0,0002	%		
	Mercury (Hg)	[7439-97-6]	0,0001	%		
	Lithium (Li)	[7439-93-2]	0,00007	%		
	Sodium (Na)	[7440-23-5]	0,00003	%		
	Lead (Pb)	[7439-92-1]	0,0035	%		