

CONOSTAN D-Series(R) D-19 Std: 30 ppm

Art. ID CON-150-300-011
Unit 100 g
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,003	%		
	Aluminium (Al)	[7429-90-5]	30	µg/g		
	Barium (Ba)	[7440-39-3]	0,003	%		
	Barium (Ba)	[7440-39-3]	30	µg/g		
	Boron (B)	[7440-42-8]	0,003	%		
	Boron (B)	[7440-42-8]	30	µg/g		
	Cadmium (Cd)	[7440-43-9]	0,003	%		
	Cadmium (Cd)	[7440-43-9]	30	µg/g		
	Chromium (Cr)	[7440-47-3]	0,003	%		
	Chromium (Cr)	[7440-47-3]	30	µg/g		
	Copper (Cu)	[7440-50-8]	0,003	%		
	Copper (Cu)	[7440-50-8]	30	µg/g		
	Iron (Fe)	[7439-89-6]	0,003	%		
	Iron (Fe)	[7439-89-6]	30	µg/g		
	Lead (Pb)	[7439-92-1]	0,003	%		
	Lead (Pb)	[7439-92-1]	30	µg/g		
	Magnesium (Mg)	[7439-95-4]	0,003	%		
	Magnesium (Mg)	[7439-95-4]	30	µg/g		
	Manganese (Mn)	[7439-96-5]	0,003	%		
	Manganese (Mn)	[7439-96-5]	30	µg/g		
	Molybdenum (Mo)	[7439-98-7]	0,003	%		
	Molybdenum (Mo)	[7439-98-7]	30	µg/g		
	Nickel (Ni)	[7440-02-0]	0,003	%		
	Nickel (Ni)	[7440-02-0]	30	µg/g		
	Silicon (Si)	[7440-21-3]	0,003	%		
	Silicon (Si)	[7440-21-3]	30	µg/g		
	Silver (Ag)	[7440-22-4]	0,003	%		
	Silver (Ag)	[7440-22-4]	30	µg/g		
	Sodium (Na)	[7440-23-5]	0,003	%		
	Sodium (Na)	[7440-23-5]	30	µg/g		
	Tin (Sn)	[7440-31-5]	0,003	%		

Tin (Sn)	[7440-31-5]	30	µg/g
Titanium (Ti)	[7440-32-6]	0,003	%
Titanium (Ti)	[7440-32-6]	30	µg/g
Vanadium (V)	[7440-62-2]	0,003	%
Vanadium (V)	[7440-62-2]	30	µg/g
Zinc (Zn)	[7440-66-6]	0,003	%
Zinc (Zn)	[7440-66-6]	30	µg/g