

(alphaR)-alpha-Azidobenzenepropanoic acid

Art. ID TRC-A919905-250MG
Unit 250 mg
Deliverydetails to be checked

Description

Application Note: (alphaR)-alpha-Azidobenzenepropanoic acid is used as a reagent in the synthesis of (R)- and (S)-alpha-azido carboxylic acids.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	(alphaR)-alpha-Azidobenzenepropanoic acid	[118460-00-7]				