

2-Amino-acetonitrile-2,2-d2

Art. ID TRC-A580097-100MG
Unit 100 mg
Deliverydetails to be checked

Description

Application Note: Isotope labelled Aminoacetonitrile (A580055), a useful synthetic intermediate. It is a reagent used to synthesize dipeptide nitriles as reversible and potent cathepsin S inhibitors. It can also used to prepare substituted cyclic ureas as possible HIV protease inhibitors.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2-Amino-acetonitrile-2, 2-d2	[1222011-50-8]				