

(R)-4-(2-Aminopropyl)phenol

Art. ID TRC-A621685-1G
Unit 1 g
Deliverydetails to be checked

Description

Application Note: (R)-4-(2-Aminopropyl)phenol is a reagent used for research. (R)-4-(2-Aminopropyl)phenol has been applied in multivariate optimization of zebrafish water tank model for investigation of biotransformation products and metabolism of sibutramine, selegiline, JWH-073 and hexareline drugs for doping control analysis.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	(R)-4-(2-Aminopropyl)ph enol	[1518-89-4]				