

## **N-Acetyl-3-methoxytyrosine**

Art. ID TRC-A459635-500MG  
Unit 500 mg  
Deliverydetails to be checked

### Description

Application Note: N-Acetyl-3-methoxytyrosine is an intermediate in the synthesis of rac 3-O-Methyl DOPA (Levodopa EP Impurity C) (M303805), which is an impurity of Levodopa (D533751). A useful marker that can be used to screen for the absence of Aromatic L-amino acid decarboxylase (AADC).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	N-Acetyl-3-methoxytyrosine	[63357-45-9]				