

ARCONIC - Aluminium base, RM, 213, wrought, Disc 64 mm x 25 mm

Art. ID ALCOA-KA-213
Unit disc
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Chromium (Cr) | [7440-47-3] | <0,0050 | % | | |
| | Copper (Cu) | [7440-50-8] | 7,0i | % | | |
| | Iron (Fe) | [7439-89-6] | 1,0i | % | | |
| | Magnesium (Mg) | [7439-95-4] | 0,05i | % | | |
| | Manganese (Mn) | [7439-96-5] | 0,30i | % | | |
| | Nickel (Ni) | [7440-02-0] | 0,20i | % | | |
| | Silicon (Si) | [7440-21-3] | 2,0i | % | | |
| | Tin (Sn) | [7440-31-5] | 0,55 | % | | |
| | Titanium (Ti) | [7440-32-6] | 0,05i | % | | |
| | Zinc (Zn) | [7440-66-6] | 1,0i | % | | |