

ARCONIC - Aluminium base, Mg/Zn, SUS, wrought, typical analysis, Disc 64 mm x 37 mm

Art. ID ALCOA-SQ-11
Unit disc
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Beryllium (Be)	[7440-41-7]	0,005	%		
	Calcium (Ca)	[7440-70-2]	0,02	%		
	Cobalt (Co)	[7440-48-4]	0,01	%		
	Chromium (Cr)	[7440-47-3]	0,25	%		
	Copper (Cu)	[7440-50-8]	0,5	%		
	Iron (Fe)	[7439-89-6]	0,2	%		
	Gallium (Ga)	[7440-55-3]	0,03	%		
	Magnesium (Mg)	[7439-95-4]	3	%		
	Manganese (Mn)	[7439-96-5]	0,4	%		
	Silicon (Si)	[7440-21-3]	0,2	%		
	Titanium (Ti)	[7440-32-6]	0,1	%		
	Zinc (Zn)	[7440-66-6]	6,6	%		