

## **0.5 mg/L C from USP Sucrose Total Organic Carbon (TOC) Standard for Anatel A643 Analyser**

Art. ID	REA-ISTOC1189
Unit	60 mL
Deliverydetails	No Dangerous Good /not restricted

### Description

Reagecon's 0.5 mg/L C from USP Sucrose Total Organic Carbon (TOC) Standard is a product of extremely high specification and purity. The product is produced from ultra pure raw materials including water that is prepared by a special proprietary process. The product is presented in bottles that are manufactured, cleaned and leached, by a special proprietary process. specifically for the product

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Total Organic Carbon (T OC)		mg/L			