

0.5 mg/L C from USP 1,4- Benzoquinone Total Organic Carbon (TOC) Standard for Anatel A643 Analyser

Art. ID REA-ISTOC1190
Unit 60 mL
Deliverydetails No Dangerous Good /not restricted

Description

Reagecon's 0.5 mg/L C from USP 1,4- Benzoquinone Total Organic Carbon (TOC) Standard is a product of extremely high specification and purity. The product is produced from ultra pure raw materials including water that is prepared by a special proprietary process. The product is presented in bottles that are manufactured, cleaned and leached, by a special proprietary process specifically for the product

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Total Organic Carbon (T OC)		mg/L			