

**Custom Solution - Selenium Se - 100 mg/l in HNO3 for ICP-MS (Se 99.999%)**

Art. ID LM24-CUS-70612  
Unit 125 mL bottle  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

**Description**

Conservation 24 months. (12 months after opening).

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Selenium (Se)     | [7782-49-2] | 100                 | mg/L |        |        |