

**Thioridazine-D3 hydrochloride solution, 100 µg/mL in methanol, certified reference material, ampule of 1 mL, Cerilliant(R) (SUPELCO)**

Art. ID SAF-T-135-1ML

Unit 1 x 1 mL

| Text/Information | Analyte/Parameter             | CAS number     | Concentration/Value | Unit | Method | Source |
|------------------|-------------------------------|----------------|---------------------|------|--------|--------|
|                  | Thioridazine-d3 Hydrochloride | [1189928-36-6] |                     |      |        |        |