

Propionyl-L-Carnitine-(N- Methyl-D3) solution, 500 µg/mL (90:10 Methanol:DMSO), certified reference material, ampule of 1 mL, Cerilliant(R) (SUPELCO)

Art. ID SAF-A-163-1ML
Unit 1 x 1 mL
Deliverydetails Dry ice shipment

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Propionyl-L-carnitine-(N-methyl-d3)	[203806-01-3]				