

Palmitoyl-L-Carnitine (N-Methyl-D3) solution, 500 µg/mL (90:10 Methanol:DMSO), certified reference material, ampule of 1 mL, Cerilliant(R) (SUPELCO)

Art. ID SAF-A-164-1ML
Unit 1 x 1 mL
Deliverydetails Dry ice shipment

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Palmitoyl-L-carnitine-(N-methyl-d3)	[202480-73-7]				