

**Palmitoyl-L-Carnitine (N-Methyl-D3) solution, 500 µg/mL (90:10 Methanol:DMSO), certified reference material, ampule of 1 mL, Cerilliant(R) (SUPELCO)**

Art. ID                      SAF-A-164-1ML  
Unit                         1 x 1 mL  
Deliverydetails            Dry ice shipment

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Palmitoyl-L-carnitine-( N-methyl-d3)	[202480-73-7]				