

Creatinine-D3 solution, 100 µg/mL ((50:50) MeOH:Water), certified reference material, ampule of 1 mL, Cerilliant(R) (SUPELCO)

Art. ID SAF-C-165-1ML
Unit 1 x 1 mL
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|---------------|---------------------|------|--------|--------|
| | Creatinine-d3 | [143827-20-7] | | | | |