

ARCONIC - Aluminium base, Cu/Mg/Zn, RM, wrought, Disc 64 mm x 25 mm

Art. ID ALCOA-ST1-7000
Unit disc
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	0,001i	%		
	Boron (B)	[7440-42-8]	<0,0005	%		
	Beryllium (Be)	[7440-41-7]	0,0005	%		
	Bismuth (Bi)	[7440-69-9]	0,001i	%		
	Cadmium (Cd)	[7440-43-9]	0,001i	%		
	Cobalt (Co)	[7440-48-4]	0,001i	%		
	Chromium (Cr)	[7440-47-3]	0,15i	%		
	Copper (Cu)	[7440-50-8]	1,5i	%		
	Iron (Fe)	[7439-89-6]	0,25	%		
	Gallium (Ga)	[7440-55-3]	0,02i	%		
	Mercury (Hg)	[7439-97-6]	0,003i	%		
	Magnesium (Mg)	[7439-95-4]	2,5i	%		
	Manganese (Mn)	[7439-96-5]	0,25	%		
	Nickel (Ni)	[7440-02-0]	0,02i	%		
	Lead (Pb)	[7439-92-1]	0,003i	%		
	Antimony (Sb)	[7440-36-0]	0,001i	%		
	Silicon (Si)	[7440-21-3]	0,15i	%		
	Tin (Sn)	[7440-31-5]	0,001i	%		
	Strontium (Sr)	[7440-24-6]	0,001i	%		
	Titanium (Ti)	[7440-32-6]	0,04i	%		
	Vanadium (V)	[7440-62-2]	0,01i	%		
	Zinc (Zn)	[7440-66-6]	6,5i	%		
	Zirconium (Zr)	[7440-67-7]	0,001i	%		