

**ARCONIC - Aluminium base, RM, 3004, wrought, Disc 64 mm x 25 mm**

Art. ID ALCOA-WB-3004  
Unit disc  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Bismuth (Bi)	[7440-69-9]	0,002i	%		
	Cadmium (Cd)	[7440-43-9]	0,003i	%		
	Chromium (Cr)	[7440-47-3]	0,02i	%		
	Copper (Cu)	[7440-50-8]	0,20i	%		
	Iron (Fe)	[7439-89-6]	0,4	%		
	Gallium (Ga)	[7440-55-3]	0,02i	%		
	Magnesium (Mg)	[7439-95-4]	0,9	%		
	Manganese (Mn)	[7439-96-5]	1,4i	%		
	Nickel (Ni)	[7440-02-0]	0,02i	%		
	Lead (Pb)	[7439-92-1]	0,02i	%		
	Silicon (Si)	[7440-21-3]	0,10i	%		
	Titanium (Ti)	[7440-32-6]	0,05i	%		
	Vanadium (V)	[7440-62-2]	0,01i	%		
	Zinc (Zn)	[7440-66-6]	0,10i	%		